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

1- Description and functions

The **DynaStar** foot includes the following main components :

- A carbon fiber blade that allows for energy storage and return, as well as inversion/eversion at forefoot.
- A light alloy frame that ensures heel strike and transfers the stresses to the blade.

These two parts are linked by an elastomer axis that provides inversion/eversion of the foot and a good transfer of weight bearing from the heel to the forefoot.

- A polyurethane cosmetic shell which transmits the mechanical stresses between the shoe and the composite structure.
 - A sock in technical fibers reduces noises and wear due to contact between the cosmetic shell and the structure.
- Thanks to these components, the foot has a comfortable and flexible behaviour.

2- Indications, maximum weight and connecting components


The **DynaStar** foot is designed exclusively for lower limb prostheses. It is appropriate for patients with low, medium or high activity level whose weight (carried load inclusive) complies with the following table.

The **DynaStar** foot is designed for a 10 mm heel height. It includes a male pyramid that allows for the use of standard connectors available on the market. Thanks to these connectors, the heel height can be easily modified from 0 to 20 mm.

Module selection guide				
	<i>Patient weight (Kg)</i>			
<i>Activity</i>	45 to 60	60 to 80	80 to 100	100 to 125
<i>Moderate</i>	2	2	3	4
<i>Normal</i>	2	3	4	5
<i>High</i>	3	4	5	N/A

3- Assembling (without the patient)

Place under the foot a wedge of the same height as the shoe heel. Fix the female pyramid connector to the foot, taking care that its axis be perpendicular to the ground (fig. 1).

4- Alignment (with the patient): Place the foot so that the load line goes through the posterior third of the foot (fig. 2)



Fig. 2

Notes:

- Increasing equinus provides a more dynamic alignment. It increases energy return at step completion.
- If the foot seems too long during transfer to the forefoot or if heel support is too flexible, you have to position the foot behind the socket.
- If heel support seems too hard, or if the foot seems too short during transfer to the forefoot, you have to position the foot ahead of the socket.

5- Finishing: To avoid rubbing noises, pull the sock up to the ankle connector or the basis of the pyramid and fix it with a fine string for example.

The malleolar clips shown in the following table ensure connection to the cosmetic foam and allow to dismantle the foot easily and rapidly.

	Foot size (cm)	22 to 25	26 to 27	28 to 29
	Product Nr	1A40013	1A40014	1A40015

Installation: Cut the cosmetic foam to the desired length, position the foam on the foot and draw the malleolus contour. Inside the foam mill a space that will receive the selected connecting part and the foot part protruding from the shell, and at the same time take care that thickness be sufficient to allow for sticking. Position the malleolar clip carefully and clip it into the cosmetic foam. Clean contact surfaces with isopropyl alcohol and stick the malleolar clip on the foam with XC090 neoprene adhesive. After drying, finish the outer shape as usually.


6- Maintenance

The cosmetic shell, the sock and the malleolar clip are the only foot components that may be worn out and may need to be replaced.

The cosmetic shell can be removed and installed again using 1A40019 shoe-horn.

Use exclusively the spare cosmesis 1A40020-D22 including the cosmesis, the sock and the malleolar clip.

Side Size in cm

7- Recycling

The different components of the foot are special wastes and must be treated accordingly.

Part for the patient

ADVICE FOR USE, MAINTENANCE AND SAFETY

Clean the cosmetic cover with a slightly damp cloth or sponge.

Use temperature : - 20°C / +60 °C.

Do not let foreign bodies or dust inside the foot.

If the foot has been dipped or sprinkled, empty the remaining liquid and dry the foot rapidly.

A use of the foot that does not adhere to your prosthetist's recommendations can damage foot components.

In case of an abnormal behaviour of your foot, or changes in your foot features, or if the foot sustained a hard shock, contact your prosthetist.

PROTEOR
6 rue de la Redoute - 21850 Saint-Apollinaire - France
Tél. 33 (0) 3 80 78 42 10 - Fax 33 (0) 3 80 78 42 15
ht.orthopedie@proteor.com - www.proteor.com

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