# Focused Solutions

Freedom Innovations is solely focused on developing worldclass lower-limb solutions in close collaboration with amputees and prosthetists. These solutions encompass technological innovations, service, training, educational resources and consultative support to help ensure user satisfaction.

## Freedom Agilix™ Technical Specifications

Traditional sizes: 22-31 cm Sandal Toe Sizes: 22-28 cm

**Build Height:** 143 mm-172 mm (5.6- 6.8 in) **Product Weight:** 19-27 oz (553-772g) **User Weight Limit:** 166 kg (365 lbs)

**Stiffness Categories:** 1-9 (low and mod impact)

Split Keel (Inversion/Eversion): Yes

Heel Height: 10 mm (3/8 in) Foot Shell CAP: Upon request

Footshell skin tones: 3 (light, medium, dark)

Attachment Type: Male pyramid

Warranty: Foot module 36 months, Foot shell 6 months

### Determine Your Impact Level

Low	Moderate	
Walking	Light jogging	
Gardening	Tennis	
Shopping	Volleyball	
Golfing	Hiking	
Biking	Dancing	
Fishing	Skating	
Camping		



Weight (lbs)	Impact Level		Weight (kg)
	Low	Moderate	
100-115	1	1	44-52
116-130	1	2	53-59
131-150	2	3	60-68
151-170	3	4	69-77
171-195	4	5	78-88
196-220	5	6	89-100
221-255	6	7	101-116
256-275	7	8	117-125
286-325	8	9	131-147
326-365	9		148-166

### PROTEOR

European Headquarters 6 rue de la Redoute 21850 St Aplollinaire France Phone: +33 (0) 80 78 42 08 Email: cs@proteor.com www.proteor.com

Connect with us



#proteorusa











**PROTEOR USA** 

Tempe, AZ 85282

info@proteorusa.com

www.proteorusa.com

Suite 101

1236 West Southern Avenue

Phone: 855.450.7300 Email:

# USA Manufacturing & Returns

425 East 400 North Gunnison, UT 84634











The Freedom Agilix™ is a multi-axial, shock absorbing foot system designed to manage loading impacts, reduce socket shear forces and improve comfort while walking on nearly any terrain. The ultra-lightweight design offers low to moderate K3 ambulators shockingly comfortable performance at a value like no other.



# Maximum Shock Value

# Target Users K<sub>3</sub>



- Users who prefer a softer, easy to roll-over, foot system
- Users who walk for extended periods for their vocation or exercise

Not recommended: High impact or K4 level users; athletes without a secondary prosthesis



